


TOSCA Laser Control System - Database Schema

- [Architecture Diagrams](#)
 - [Figure 1: TOSCA Database Schema ERD](#)
 - [Figure 2: TOSCA Data Architecture](#)
- [Overview](#)
- [Entity Relationship Diagram](#)
- [Table Definitions](#)
 - [1. tech_users](#)
 - [2. subjects](#)
 - [3. protocols](#)
 - [4. sessions](#)
 - [5. treatment_events](#)
 - [6. session_recordings](#)
 - [7. safety_log](#)
 - [8. calibrations](#)
 - [9. system_config](#)
- [Views \(Convenience Queries\)](#)
 - [View: Active Sessions](#)
 - [View: Subject Session History](#)
 - [View: Safety Events Summary](#)
- [Common Queries](#)
 - [1. Find Subject by Code](#)
 - [2. Get All Sessions for Subject](#)
 - [3. Get Treatment Events for Session](#)
 - [4. Get Active Calibration for Device](#)
 - [5. Safety Events in Last 24 Hours](#)
 - [6. Protocol Usage Statistics](#)
- [Database Initialization Script](#)
- [Migration Strategy](#)
- [Backup Strategy](#)
- [Performance Considerations](#)
- [Data Integrity](#)

Architecture Diagrams

Figure 1: TOSCA Database Schema ERD



TOSCA Database Schema ERD

Figure 2: TOSCA Data Architecture


```

        FOREIGN KEY (created_by_tech_id) REFERENCES tech_users(tech_id)
    );

-- Indexes
CREATE UNIQUE INDEX idx_subject_code ON subjects(subject_code);
CREATE INDEX idx_subject_active ON subjects(is_active);
CREATE INDEX idx_subject_created_date ON subjects(created_date);

```

Subject Code Format: P-YYYY-NNNN (e.g., P-2025-0001)

Sample Data:

```

INSERT INTO subjects (subject_code, created_by_tech_id, notes) VALUES
('P-2025-0001', 1, 'Initial consultation');

```

3. protocols

Saved treatment protocols (templates) that can be reused.

```

CREATE TABLE protocols (
    protocol_id INTEGER PRIMARY KEY AUTOINCREMENT,
    protocol_name TEXT UNIQUE NOT NULL,
    description TEXT,
    protocol_type TEXT NOT NULL, -- 'constant', 'linear_ramp', 'multi_step', 'custom'
    protocol_data TEXT NOT NULL, -- JSON: [{step: 1, duration: 60, power_start: 5.0, ...}, ...]
    created_date DATETIME DEFAULT CURRENT_TIMESTAMP,
    created_by_tech_id INTEGER,
    last_modified_date DATETIME,
    last_used_date DATETIME,
    usage_count INTEGER DEFAULT 0,
    is_active BOOLEAN DEFAULT 1,

    FOREIGN KEY (created_by_tech_id) REFERENCES tech_users(tech_id)
);

-- Indexes
CREATE UNIQUE INDEX idx_protocol_name ON protocols(protocol_name);
CREATE INDEX idx_protocol_type ON protocols(protocol_type);
CREATE INDEX idx_protocol_active ON protocols(is_active);

```

Protocol Data JSON Format:

```

{
  "version": "1.0",
  "steps": [
    {
      "step_number": 1,
      "duration_seconds": 60,
      "power_start_watts": 5.0,
      "power_end_watts": 5.0,
      "ring_size_mm": 3.0,
      "ramp_type": "constant"
    },
    {
      "step_number": 2,
      "duration_seconds": 30,
      "power_start_watts": 5.0,
      "power_end_watts": 8.0,
      "ring_size_mm": 3.5,
      "ramp_type": "linear"
    }
  ],
  "safety_limits": {
    "max_power_watts": 10.0,
    "max_duration_seconds": 300
  }
}

```

4. sessions

Treatment sessions - each subject visit creates one session.

```

CREATE TABLE sessions (
    session_id INTEGER PRIMARY KEY AUTOINCREMENT,
    subject_id INTEGER NOT NULL,
    tech_id INTEGER NOT NULL,
    protocol_id INTEGER, -- NULL if custom/ad-hoc treatment

    start_time DATETIME NOT NULL,
    end_time DATETIME,
    duration_seconds INTEGER, -- Calculated: end_time - start_time

    status TEXT NOT NULL, -- 'in_progress', 'completed', 'aborted', 'paused'
    abort_reason TEXT, -- Filled if status = 'aborted'

    treatment_site TEXT, -- Anatomical location
    treatment_notes TEXT, -- Pre-treatment notes

    protocol_name TEXT, -- Snapshot of protocol name (in case protocol deleted)
    protocol_data_snapshot TEXT, -- JSON snapshot of protocol used

    session_folder_path TEXT, -- Path to session data folder
    video_path TEXT, -- Path to recorded video

    -- Summary statistics (calculated at end)
    total_laser_on_time_seconds REAL,
    avg_power_watts REAL,
    max_power_watts REAL,
    total_energy_joules REAL,

    post_treatment_notes TEXT, -- Post-treatment observations
    created_date DATETIME DEFAULT CURRENT_TIMESTAMP,

    FOREIGN KEY (subject_id) REFERENCES subjects(subject_id),
    FOREIGN KEY (tech_id) REFERENCES tech_users(tech_id),
    FOREIGN KEY (protocol_id) REFERENCES protocols(protocol_id)
);

```

```
-- Indexes
CREATE INDEX idx_session_subject ON sessions(subject_id);
CREATE INDEX idx_session_tech ON sessions(tech_id);
CREATE INDEX idx_session_start_time ON sessions(start_time);
CREATE INDEX idx_session_status ON sessions(status);
CREATE INDEX idx_session_date_status ON sessions(start_time, status);
```

5. treatment_events

High-frequency event log - every significant action during treatment.

Note: For very high-frequency data (>10Hz), consider writing to JSON file first, then importing to DB post-session.

```
CREATE TABLE treatment_events (
  event_id INTEGER PRIMARY KEY AUTOINCREMENT,
  session_id INTEGER NOT NULL,
  timestamp DATETIME NOT NULL,
  time_offset_ms INTEGER NOT NULL, -- Milliseconds from session start

  event_type TEXT NOT NULL,
  -- Event types:
  -- 'laser_on', 'laser_off'
  -- 'power_set', 'power_change'
  -- 'ring_size_set', 'ring_size_change'
  -- 'protocol_step_start', 'protocol_step_end'
  -- 'pause', 'resume'
  -- 'image_captured'
  -- 'safety_warning', 'safety_trigger'
  -- 'user_action'

  -- Hardware state at time of event
  laser_power_commanded_watts REAL,
  laser_power_actual_watts REAL, -- From photodiode
  ring_size_mm REAL,
  actuator_position_um REAL,
  photodiode_voltage_v REAL,

  -- Safety status
  footpedal_state BOOLEAN,
  smoothing_device_state BOOLEAN,
  interlock_status TEXT, -- JSON: {"footpedal": true, "smoothing": true, ...}

  -- Image/video reference
  image_path TEXT,
  video_timestamp_ms INTEGER,

  -- Additional data (JSON for flexibility)
  event_details TEXT, -- JSON: free-form additional data

  FOREIGN KEY (session_id) REFERENCES sessions(session_id)
);

-- Indexes
CREATE INDEX idx_event_session ON treatment_events(session_id);
CREATE INDEX idx_event_timestamp ON treatment_events(timestamp);
CREATE INDEX idx_event_type ON treatment_events(event_type);
CREATE INDEX idx_event_session_time ON treatment_events(session_id, time_offset_ms);
```

Event Details JSON Examples:

```
// laser_on event
{
  "source": "footpedal_depression",
  "protocol_step": 1
}

// safety_trigger event
{
  "trigger_source": "photodiode_out_of_range",
  "expected_power": 5.0,
  "measured_power": 2.3,
  "action_taken": "laser_shutdown"
}

// user_action event
{
  "action": "protocol_power_override",
  "old_value": 5.0,
  "new_value": 6.5,
  "reason": "subject tolerance"
}
```

6. session_recordings

Metadata for recorded video and image files.

```
CREATE TABLE session_recordings (
  recording_id INTEGER PRIMARY KEY AUTOINCREMENT,
  session_id INTEGER NOT NULL,
  recording_type TEXT NOT NULL, -- 'video', 'snapshot', 'focus_stack'
  file_path TEXT NOT NULL,
  file_size_bytes INTEGER,
  duration_seconds REAL, -- For videos
  timestamp DATETIME NOT NULL,
  frame_count INTEGER, -- For videos
  resolution TEXT, -- e.g., '1920x1080'
  notes TEXT,

  FOREIGN KEY (session_id) REFERENCES sessions(session_id)
);

-- Indexes
CREATE INDEX idx_recording_session ON session_recordings(session_id);
CREATE INDEX idx_recording_type ON session_recordings(recording_type);
```

7. safety_log

Comprehensive log of all safety-related events (separate from treatment_events for clarity).

```
CREATE TABLE safety_log (
  log_id INTEGER PRIMARY KEY AUTOINCREMENT,
  timestamp DATETIME NOT NULL,
  session_id INTEGER, -- NULL if not during active session
  tech_id INTEGER,

  event_type TEXT NOT NULL,
  -- Event types:
  -- 'e_stop_pressed', 'e_stop_released'
  -- 'footpedal_released', 'footpedal_timeout'
  -- 'smoothing_device_fault', 'smoothing_device_recovery'
  -- 'photodiode_out_of_range', 'photodiode_fault'
  -- 'interlock_failure', 'interlock_recovery'
  -- 'power_limit_exceeded'
  -- 'system_watchdog_timeout'

  severity TEXT NOT NULL, -- 'info', 'warning', 'critical', 'emergency'
  description TEXT NOT NULL,

  -- Hardware state at time of event
  system_state TEXT, -- 'idle', 'ready', 'armed', 'treating', 'fault'
  laser_state TEXT,
  footpedal_state BOOLEAN,
  smoothing_device_state BOOLEAN,
  photodiode_voltage_v REAL,

  action_taken TEXT, -- What the system did in response
  details TEXT, -- JSON: additional data

  FOREIGN KEY (session_id) REFERENCES sessions(session_id),
  FOREIGN KEY (tech_id) REFERENCES tech_users(tech_id)
);

-- Indexes
CREATE INDEX idx_safety_timestamp ON safety_log(timestamp);
CREATE INDEX idx_safety_session ON safety_log(session_id);
CREATE INDEX idx_safety_severity ON safety_log(severity);
CREATE INDEX idx_safety_type ON safety_log(event_type);
```

8. calibrations

Device calibration data and history.

```
CREATE TABLE calibrations (
  calibration_id INTEGER PRIMARY KEY AUTOINCREMENT,
  calibration_date DATETIME NOT NULL,
  tech_id INTEGER NOT NULL,

  device_type TEXT NOT NULL, -- 'laser', 'actuator', 'photodiode', 'camera'
  device_serial_number TEXT,

  calibration_type TEXT NOT NULL,
  -- Calibration types:
  -- 'actuator_position_to_ring_size' -- Maps actuator µm to ring mm
  -- 'photodiode_voltage_to_power' -- Maps voltage to laser watts
  -- 'laser_power_verification' -- Verifies commanded vs actual power
  -- 'camera_focus_calibration' -- Focus score reference
  -- 'ring_detection_calibration' -- Circle detection tuning

  calibration_data TEXT NOT NULL, -- JSON: calibration parameters
  validation_data TEXT, -- JSON: validation results

  is_active BOOLEAN DEFAULT 1, -- Most recent calibration is active
  notes TEXT,

  FOREIGN KEY (tech_id) REFERENCES tech_users(tech_id)
);

-- Indexes
CREATE INDEX idx_calib_device ON calibrations(device_type);
CREATE INDEX idx_calib_date ON calibrations(calibration_date);
CREATE INDEX idx_calib_active ON calibrations(is_active);
```

Calibration Data JSON Examples:

```
// actuator_position_to_ring_size
{
  "calibration_points": [
    {"actuator_position_um": 0, "ring_size_mm": 2.0},
    {"actuator_position_um": 1000, "ring_size_mm": 3.0},
    {"actuator_position_um": 2000, "ring_size_mm": 4.0}
  ],
  "interpolation_method": "linear",
  "fit_equation": "ring_size = 2.0 + (position / 1000)"
}

// photodiode_voltage_to_power
{
  "calibration_points": [
    {"voltage_v": 0.5, "power_watts": 1.0},
    {"voltage_v": 2.5, "power_watts": 5.0},
    {"voltage_v": 5.0, "power_watts": 10.0}
  ],
  "fit_equation": "power = voltage * 2.0",
  "r_squared": 0.998
}
```

9. system_config

Application configuration and settings (key-value store).

```
CREATE TABLE system_config (
  config_key TEXT PRIMARY KEY,
  config_value TEXT NOT NULL,
  config_type TEXT NOT NULL, -- 'string', 'integer', 'float', 'boolean', 'json'
```

```

description TEXT,
last_modified DATETIME DEFAULT CURRENT_TIMESTAMP,
modified_by_tech_id INTEGER,

FOREIGN KEY (modified_by_tech_id) REFERENCES tech_users(tech_id)
);

-- Sample configuration entries
INSERT INTO system_config (config_key, config_value, config_type, description) VALUES
('max_laser_power_watts', '10.0', 'float', 'Hard limit for laser power'),
('treatment_loop_frequency_hz', '10', 'integer', 'Control loop update rate'),
('photodiode_warning_threshold_percent', '15', 'float', 'Deviation % before warning'),
('photodiode_fault_threshold_percent', '30', 'float', 'Deviation % before fault'),
('video_resolution', '1920x1080', 'string', 'Recording resolution'),
('video_fps', '30', 'integer', 'Recording frame rate'),
('session_data_path', 'data/sessions', 'string', 'Base path for session data');

```

Views (Convenience Queries)

View: Active Sessions

```

CREATE VIEW active_sessions AS
SELECT
    s.session_id,
    s.subject_id,
    p.subject_code,
    s.tech_id,
    t.full_name as tech_name,
    s.start_time,
    s.status,
    s.protocol_name,
    strftime('%s', 'now') - strftime('%s', s.start_time) as duration_seconds
FROM sessions s
JOIN subjects p ON s.subject_id = p.subject_id
JOIN tech_users t ON s.tech_id = t.tech_id
WHERE s.status IN ('in_progress', 'paused');

```

View: Subject Session History

```

CREATE VIEW subject_session_history AS
SELECT
    p.subject_id,
    p.subject_code,
    s.session_id,
    s.start_time,
    s.end_time,
    s.status,
    s.protocol_name,
    s.total_laser_on_time_seconds,
    s.avg_power_watts,
    s.total_energy_joules,
    t.full_name as technician
FROM subjects p
LEFT JOIN sessions s ON p.subject_id = s.subject_id
LEFT JOIN tech_users t ON s.tech_id = t.tech_id
ORDER BY p.subject_id, s.start_time DESC;

```

View: Safety Events Summary

```

CREATE VIEW safety_events_summary AS
SELECT
    DATE(timestamp) as event_date,
    severity,
    event_type,
    COUNT(*) as event_count
FROM safety_log
GROUP BY DATE(timestamp), severity, event_type
ORDER BY event_date DESC, severity;

```

Common Queries

1. Find Subject by Code

```

SELECT * FROM subjects
WHERE subject_code = 'P-2025-0001';

```

2. Get All Sessions for Subject

```

SELECT
    s.session_id,
    s.start_time,
    s.end_time,
    s.status,
    s.protocol_name,
    s.total_energy_joules,
    t.full_name as technician
FROM sessions s
JOIN tech_users t ON s.tech_id = t.tech_id
WHERE s.subject_id = ?
ORDER BY s.start_time DESC;

```

3. Get Treatment Events for Session

```

SELECT
    timestamp,
    event_type,
    laser_power_commanded_watts,
    laser_power_actual_watts,
    ring_size_mm,
    photodiode_voltage_v
FROM treatment_events
WHERE session_id = ?
ORDER BY time_offset_ms;

```

4. Get Active Calibration for Device

```
SELECT calibration_data
FROM calibrations
WHERE device_type = 'photodiode'
      AND is_active = 1
ORDER BY calibration_date DESC
LIMIT 1;
```

5. Safety Events in Last 24 Hours

```
SELECT
    timestamp,
    event_type,
    severity,
    description,
    session_id
FROM safety_log
WHERE timestamp >= datetime('now', '-1 day')
ORDER BY timestamp DESC;
```

6. Protocol Usage Statistics

```
SELECT
    protocol_name,
    usage_count,
    last_used_date,
    (SELECT COUNT(*) FROM sessions WHERE protocol_id = p.protocol_id) as actual_usage
FROM protocols p
WHERE is_active = 1
ORDER BY usage_count DESC;
```

Database Initialization Script

```
-- Enable foreign key constraints
PRAGMA foreign_keys = ON;

-- Set journal mode for better concurrency
PRAGMA journal_mode = WAL;

-- Create all tables (in dependency order)
-- ... (include all CREATE TABLE statements above) ...

-- Create views
-- ... (include all CREATE VIEW statements above) ...

-- Insert default tech user
INSERT INTO tech_users (username, full_name, role) VALUES
    ('admin', 'System Administrator', 'admin');

-- Insert default system config
-- ... (include system_config INSERTs) ...

-- Log database creation
INSERT INTO safety_log (timestamp, event_type, severity, description, system_state, action_taken)
VALUES (CURRENT_TIMESTAMP, 'database_initialized', 'info', 'Database schema created', 'idle', 'none');
```

Migration Strategy

Use Alembic for schema migrations:

```
# alembic/versions/001_initial_schema.py
def upgrade():
    # Create tables
    pass

def downgrade():
    # Rollback
    pass
```

Versioning convention: YYYYMMDD_NNN_description.py

Backup Strategy

- 1. **Automatic backups:** Before each session, backup database
- 2. **Backup location:** data/backups/laser_control_YYYYMMDD_HHMMSS.db
- 3. **Retention:** Keep 30 days of daily backups
- 4. **Session data:** Backed up separately (video files are large)

Performance Considerations

- 1. **Indexing:** All foreign keys and frequently queried columns indexed
- 2. **Partitioning:** Consider separate tables for archived sessions (>1 year old)
- 3. **Write optimization:** Use transactions for bulk inserts (treatment_events)
- 4. **Query optimization:** Use EXPLAIN QUERY PLAN to optimize slow queries

Data Integrity

- 1. **Foreign key constraints:** Enabled (PRAGMA foreign_keys = ON)
- 2. **Write-ahead logging:** Enabled (PRAGMA journal_mode = WAL)
- 3. **Validation:** Enforce in application layer before DB write
- 4. **Immutability:** treatment_events and safety_log are append-only