

# TOSCA Safety Shutdown Policy

- [Architecture Diagrams](#)
- [Selective Shutdown Principle](#)
  - [Core Policy](#)
- [Shutdown Hierarchy](#)
- [System-by-System Shutdown Rules](#)
  - [Treatment Laser \(SHUTDOWN\)](#)
  - [Camera System \(MAINTAIN\)](#)
  - [Actuator System \(MAINTAIN\)](#)
  - [Aiming Laser \(MAINTAIN\)](#)
  - [Laser Spot Smoothing Module \(CONDITIONAL\)](#)
  - [GPIO Monitoring \(MAINTAIN\)](#)
- [Implementation in Code](#)
  - [Safety Manager Emergency Stop](#)
  - [Watchdog Failure Handler](#)
  - [Protocol Engine Safety Check](#)
- [Rationale for Selective Shutdown](#)
  - [Patient Safety Benefits](#)
  - [Diagnostic Benefits](#)
  - [Operational Benefits](#)
  - [Safety Validation](#)
- [Recovery Procedures](#)
  - [After Selective Shutdown](#)
  - [System State Transitions](#)
- [Testing Requirements](#)
  - [Selective Shutdown Tests](#)
- [Related Documents](#)
- [Revision History](#)

## Architecture Diagrams

**Document Version:** 1.1 **Last Updated:** 2025-10-26 **Status:** Active Policy **Priority:** CRITICAL

---

## Selective Shutdown Principle

### Core Policy

When safety interlocks fail or unsafe conditions are detected, **only the treatment laser** is immediately disabled. All other systems remain operational to support safe assessment and recovery.

This policy supersedes any previous documentation indicating full system shutdown.

---

## Shutdown Hierarchy

1. **\*\*Safety Interlock Failure Detected\*\***
2. **\*\*IMMEDIATE\*\*** - Treatment Laser Shutdown
3. **\*\*• Laser output OFF\*\***
4. **\*\*• Laser power 0 mW\*\***
5. **\*\*• Laser driver DISABLED\*\***
6. **\*\*• Response time\*\*** - <100ms
7. **\*\*MAINTAINED\*\*** - Support Systems
8. **\*\*[DONE] Camera (visual monitoring)\*\***

9. **\*\*[DONE]** Actuator (position control)\*\*
  10. **\*\*[DONE]** Aiming Laser (low-power alignment)\*\*
  11. **\*\*[DONE]** Laser Spot Smoothing Module (if not fault src)\*\*
  12. **\*\*[DONE]** GPIO monitoring\*\*
  13. **\*\*[DONE]** Event logging\*\*
  14. **\*\*[DONE]** UI responsiveness\*\*
- 

## System-by-System Shutdown Rules

### Treatment Laser (SHUTDOWN)

**Action:** Immediate disable **Methods:**

```
laser_controller.set_output(False) # Disable output
laser_controller.set_current(0.0)  # Zero power
```

**Rationale:** - Primary hazard source - Therapeutic power levels (up to 2000 mW) - Immediate shutdown prevents tissue damage

**Recovery:** Requires manual safety clearance and system reset

---

### Camera System (MAINTAIN)

**Action:** Continue streaming **Rationale:** - **Essential for post-fault assessment:** Operator needs to see treatment area - **Visual verification:** Confirms laser is actually off - **Diagnostic value:** May show what caused the fault - **No safety hazard:** Passive imaging only

**Example Use Cases:** - Verify laser spot has disappeared - Check for unexpected tissue reactions - Document fault conditions - Guide repositioning if needed

---

### Actuator System (MAINTAIN)

**Action:** Remain controllable **Rationale:** - **Safe repositioning:** May need to move away from treatment site - **Controlled shutdown:** Allows precise parking position - **Diagnostic access:** Can move to verify mechanical issues - **Limited hazard:** Motion at controlled speeds, low force

**Restrictions After Fault:** - Manual control only (no automatic protocols) - Reduced speed limits (safety mode) - Position bounds enforced - Continuous operator supervision required

**Example Use Cases:** - Move away from tissue after emergency stop - Return to home position - Clear workspace for fault investigation - Position for camera inspection

---

### Aiming Laser (MAINTAIN)

**Action:** Keep available **Rationale:** - **Low power:** Typically <5 mW (Class 2, inherently safe) - **Alignment verification:** Can check optical path alignment - **Visual reference:** Shows where treatment would occur - **Diagnostic tool:** Helps identify optical misalignment faults

**Safety Notes:** - Class 2 laser (eye-safe with blink reflex) - Different controller from treatment laser - Independent power supply - Can be manually disabled if desired

**Example Use Cases:** - Verify optical alignment after fault - Check for obstruction in beam path - Confirm targeting accuracy - Diagnose pointing stability issues

---

### Laser Spot Smoothing Module (CONDITIONAL)

**Action:** Depend on fault source

**If Laser Spot Smoothing Module IS the fault source:** - Motor → OFF - Safety interlock recognizes motor failure - Treatment laser disabled (per policy)

**If Laser Spot Smoothing Module is NOT the fault source:** - Motor continues running if already active - Can be manually controlled - Helps maintain component temperature stability - May aid in diagnosis

**Rationale:** - Device failure triggers selective shutdown (laser only) - Device can continue if not implicated in fault - Thermal stability may aid recovery - No direct safety hazard from motor operation

---

## GPIO Monitoring (MAINTAIN)

**Action:** Continue all monitoring functions **Rationale:** - Heartbeat monitoring must continue (watchdog still active) - photodiode laser pickoff measurement feedback confirms laser is off - Vibration sensor shows motor state - Essential for diagnosing root cause

**Monitored Signals:** - photodiode laser pickoff measurement voltage → Should drop to <0.1V after laser shutdown - Smoothing device vibration → Continues if motor running - Watchdog heartbeat → Must continue for system health - Safety interlock states → Track for diagnostics

---

## Implementation in Code

### Safety Manager Emergency Stop

```
class SafetyManager:
    def emergency_stop(self, reason: str):
        """
        Selective emergency shutdown - TREATMENT LASER ONLY.

        Maintains operational: Camera, Actuator, Aiming Laser, GPIO monitoring
        """
        logger.critical(f"EMERGENCY STOP: {reason}")

        # 1. IMMEDIATE: Disable treatment laser only
        if self.laser_controller:
            self.laser_controller.set_output(False)
            self.laser_controller.set_current(0.0)
            logger.info("Treatment laser disabled (selective shutdown)")

        # 2. MAINTAIN: Other systems remain operational
        # - Camera continues streaming
        # - Actuator remains controllable
        # - Aiming laser stays available
        # - GPIO monitoring continues
        # - Event logging active

        # 3. Update system state
        self.system_state = SystemState.FAULT
        self.laser_enabled = False

        # 4. Log event
        self.event_logger.log_event(
            event_type=EventType.EMERGENCY_STOP,
            description=f"Emergency stop: {reason} (selective shutdown - laser only)",
            severity=EventSeverity.EMERGENCY
        )

        # 5. Notify UI
        self.safety_event.emit("emergency_stop", reason)
```

### Watchdog Failure Handler

```
class SafetyWatchdog:
```

```

def _handle_critical_failure(self):
    """
    Watchdog failure → Selective shutdown (treatment laser only).

    IMPORTANT: Does not shutdown entire system.
    """
    logger.critical("Watchdog critical failure - disabling TREATMENT LASER only")

    # 1. Disable treatment laser (selective shutdown)
    if self.laser_controller:
        self.laser_controller.set_output(False)
        self.laser_controller.set_current(0.0)

    # 2. Other systems maintained:
    # - Camera: Continues monitoring
    # - Actuator: Remains controllable
    # - Aiming laser: Available
    # - GPIO: Heartbeat stops but monitoring continues

    # 3. Stop watchdog heartbeat (connection lost)
    self.stop()

    # 4. Emit critical signal
    self.watchdog_timeout_detected.emit()

    # 5. Log event
    self.event_logger.log_event(
        event_type=EventType.WATCHDOG_FAILURE,
        description="Watchdog failure - treatment laser disabled (selective shutdown)",
        severity=EventSeverity.EMERGENCY
    )

```

## Protocol Engine Safety Check

```

class ProtocolEngine:
    def _on_safety_failure_during_protocol(self):
        """
        Real-time safety failure during protocol execution.

        Selective shutdown: Stops laser, maintains monitoring.
        """
        logger.critical("Safety interlock failure during protocol execution")

        # 1. Stop protocol execution
        self.stop() # Stops protocol, sets state to STOPPED

        # 2. Disable treatment laser (selective shutdown)
        if self.laser:
            self.laser.set_output(False)
            self.laser.set_current(0.0)

        # 3. Maintain support systems:
        # - Camera continues (operator can see what happened)
        # - Actuator controllable (can reposition if needed)
        # - Event log captures full context

        # 4. Log event
        self.event_logger.log_event(
            event_type=EventType.PROTOCOL_SAFETY_STOP,
            description="Protocol stopped due to safety interlock failure (selective shutdown)",
            severity=EventSeverity.CRITICAL
        )

```

---

## Rationale for Selective Shutdown

### Patient Safety Benefits

1. **Visual Confirmation:** Operator can immediately see treatment area after fault
2. **Repositioning:** Can move away from target safely if needed

3. **Assessment:** Camera shows tissue state and reaction
4. **Documentation:** Visual record of fault conditions

## Diagnostic Benefits

1. **Root Cause Analysis:** Systems remain available for testing
2. **State Preservation:** Actuator position, camera view maintained
3. **Signal Monitoring:** GPIO continues reading sensor states
4. **Reproducibility:** Can attempt to reproduce fault with laser disabled

## Operational Benefits

1. **Reduced Downtime:** Don't need to restart entire system
2. **Faster Recovery:** Selective re-enable after fault clearance
3. **Better UX:** Operator retains control and visibility
4. **Training:** Safer environment for operator training

## Safety Validation

**Key Principle:** Treatment laser is the **only** hazard requiring immediate shutdown.

**System Risk Analysis:** - **Camera:** No harm possible (passive imaging) - **Actuator:** Low force, controlled motion, operator supervised - **Aiming Laser:** Class 2, eye-safe, inherently low risk - **Laser Spot Smoothing Module:** Mechanical device, no direct patient contact - **Treatment Laser:** High power (2000 mW), tissue-damaging, **must** be disabled

**Conclusion:** Selective shutdown optimizes both safety and functionality.

---

## Recovery Procedures

### After Selective Shutdown

**Immediate Actions:** 1. [DONE] Verify treatment laser OFF (check photodiode laser pickoff measurement < 0.1V) 2. [DONE] Assess patient/tissue status visually (camera) 3. [DONE] Document fault conditions in event log 4. [DONE] Determine root cause

**Before Re-enabling Treatment Laser:** 1. Clear the fault condition 2. Verify all safety interlocks pass 3. Confirm support systems operational 4. Reset safety manager state 5. Supervisor authorization required 6. Document recovery in event log

### System State Transitions

FAULT (Laser Disabled)

Camera: Operational (monitoring)  
Actuator: Operational (controllable)  
Aiming Laser: Operational (available)  
GPIO: Operational (monitoring)  
Event Log: Operational (recording)



[Fault Investigation & Clearance]



[Supervisor Authorization]



[Safety Checks Pass]



READY (All Systems Operational, Laser Can Be Enabled)

---

# Testing Requirements

## Selective Shutdown Tests

**Test 1:** Emergency Stop Button - Press E-stop during treatment - **Verify:** Laser OFF, camera ON, actuator controllable - **Verify:** photodiode laser pickoff measurement reads <0.1V - **Verify:** UI shows fault state but remains responsive

**Test 2:** Watchdog Timeout - Simulate GUI freeze - **Verify:** Treatment laser OFF after 1 second - **Verify:** Camera still streaming - **Verify:** Event log records fault

**Test 3:** Safety Interlock Failure (Laser Spot Smoothing Module) - Stop laser spot smoothing module during treatment - **Verify:** Laser OFF immediately - **Verify:** Motor stops, but camera/actuator available - **Verify:** Operator can reposition

**Test 4:** Selective Recovery - After fault, clear root cause - **Verify:** Can re-enable laser with supervisor auth - **Verify:** All interlocks must pass before re-enable - **Verify:** Event log documents full recovery sequence

---

## Related Documents

- docs/CODE\_REVIEW\_2025-10-26.md - Code review findings
  - docs/architecture/03\_safety\_system.md - Overall safety architecture
  - docs/architecture/06\_safety\_watchdog.md - Watchdog implementation
  - presubmit/reviews/plans/IMPLEMENTATION\_PLAN\_WATCHDOG.md - Original plan
- 

## Revision History

Version	Date	Change	Author
1.0	2025-10-15	Initial safety architecture	TOSCA Team
1.1	2025-10-26	Selective shutdown policy defined	Code Review

---

**Policy Owner:** Safety Engineer **Approval Required:** Project Lead, Safety Committee **Review Frequency:** Before each development phase **Next Review:** After Phase 1 implementation