

Calibration Procedures

- [Architecture Diagrams](#)
- [Table of Contents](#)
- [Overview](#)
 - [Purpose](#)
 - [Calibration Objectives](#)
- [Calibration Requirements](#)
 - [Calibration Hierarchy](#)
 - [Components Requiring Calibration](#)
- [Photodiode Calibration](#)
 - [Purpose](#)
 - [Equipment Required](#)
 - [Procedure](#)
- [Actuator Calibration](#)
 - [Purpose](#)
 - [Equipment Required](#)
 - [Procedure](#)
- [Camera Calibration](#)
 - [Purpose](#)
 - [Equipment Required](#)
 - [Procedure](#)
- [Power Measurement Verification](#)
 - [Purpose](#)
 - [Procedure \(Weekly\)](#)
- [Calibration Schedule](#)
 - [Routine Calibration](#)
 - [Calibration Records](#)
- [Documentation & Audit Trail](#)
 - [FDA Requirements](#)
 - [Calibration Certificate Example](#)
- [References](#)
 - [Standards](#)

Architecture Diagrams

Document Version: 1.0 **Last Updated:** 2025-10-26 **Status:** WARNING: **NOT IMPLEMENTED** - Planning Only (Phase 6+) **Priority:** CRITICAL - Required for treatment accuracy and FDA compliance

WARNING: WARNING - No Calibration Procedures Implemented:

Current Status: ALL calibration procedures documented here are **PLANNED ONLY**.

NO calibration workflows are implemented in the current codebase: - [FAILED] Photodiode calibration - NOT implemented - [FAILED] Actuator calibration - NOT implemented - [FAILED] Camera pixel calibration - NOT implemented - [FAILED] Power verification - NOT implemented - [FAILED] Calibration database storage - NOT implemented

Current Operation: System operates with factory defaults (no field calibration)

Implementation Target: Phase 6 (before clinical testing)

This document serves as: - Design specification for future calibration workflows - FDA submission documentation (planned procedures) - NIST traceability requirements documentation

Table of Contents

- [1. Overview](#)
 - [2. Calibration Requirements](#)
 - [3. Photodiode Calibration](#)
 - [4. Actuator Calibration](#)
 - [5. Camera Calibration](#)
 - [6. Power Measurement Verification](#)
 - [7. Calibration Schedule](#)
 - [8. Documentation & Audit Trail](#)
-

Overview

Purpose

TOSCA requires regular calibration of all measurement systems to ensure treatment accuracy and patient safety. This document defines calibration procedures, schedules, and documentation requirements.

Regulatory Context: - **FDA 21 CFR Part 820:** Quality System Regulation (Design Controls) - **ISO 13485:** Medical devices — Quality management systems - **IEC 60601-1:** Medical electrical equipment safety - **ISO 10012:** Measurement management systems

Calibration Objectives

- Accuracy:** Ensure all measurements are within specification
 - Traceability:** Link calibration to NIST standards
 - Safety:** Verify safety systems function correctly
 - Compliance:** Meet FDA requirements for calibration
 - Audit Trail:** Document all calibration activities
-

Calibration Requirements

Calibration Hierarchy

NIST Traceable Standards
(National Institute of Standards and Technology)

Calibration Certificate



- Reference Equipment (Certified)
- Laser power meter (certified annually)
 - Micrometer (certified annually)
 - Oscilloscope (certified every 2 years)

Used to calibrate TOSCA



- TOSCA Components
- Photodiode (power measurement)
 - Actuator (position accuracy)
 - Camera (focus measurement)
 - GPIO timing (watchdog intervals)

Components Requiring Calibration

Component	Parameter	Tolerance	Frequency	Method
Photodiode	Power measurement	±5%	Monthly	Certified power meter
Actuator	Position accuracy	±10 µm	Quarterly	Micrometer + reference block
Camera	Pixel calibration	±1 pixel	Annually	Calibration target
Laser	Output power	±10%	Weekly	Photodiode feedback loop
Watchdog	Timeout interval	±50 ms	Annually	Oscilloscope

Photodiode Calibration

Purpose

Verify photodiode accurately measures laser power (5W range).

Safety Critical: YES - Incorrect power measurement could exceed safety limits

Equipment Required

1. **Certified laser power meter** (Ophir or Thorlabs, NIST traceable)
2. **Laser controller** (Arroyo 1064nm laser)
3. **TOSCA system** (photodiode under test)

Procedure

Step 1: Setup

1. Connect certified power meter to laser output
2. Connect TOSCA photodiode to secondary measurement port
3. Allow laser to warm up for 10 minutes
4. Zero both power meters with laser off

Step 2: Calibration Points

Test at 5 power levels across range:

Target Power (W)	Certified Meter (W)	TOSCA Photodiode (V)	Acceptable Range
0.0	0.00 ± 0.01	0.00 ± 0.05	-0.05 to +0.05 V
2.5	2.50 ± 0.05	TBD	±5% of certified
5.0	5.00 ± 0.05	TBD	±5% of certified
7.5	7.50 ± 0.05	TBD	±5% of certified
10.0	10.00 ± 0.10	TBD	±5% of certified

Step 3: Calibration Curve

```
1. ## Calculate calibration curve (voltage watts)**
2. **def calibrate_photodiode(calibration_points** - List[Tuple[float, float]]) -> Tuple[float, float]:
3. ****
4. **Calculate calibration curve from voltage/power pairs.**
5. **Args** -
6. **calibration_points** - List of (voltage_v, power_w) tuples
7. **Returns** -
8. **(slope, intercept) for power = slope * voltage + intercept**
9. ****
```

```

10. import numpy as np
11. voltages = np.array([v for v, p in calibration_points])
12. powers = np.array([p for v, p in calibration_points])
13. # Linear regression
14. slope, intercept = np.polyfit(voltages, powers, deg=1)
15. return (slope, intercept)
16. # Example usage
17. cal_points = [
18.     (0.0, 0.0),    # 0V 0W
19.     (0.5, 2.5),    # 0.5V 2.5W
20.     (1.0, 5.0),    # 1.0V 5.0W
21.     (1.5, 7.5),    # 1.5V 7.5W
22.     (2.0, 10.0),   # 2.0V 10.0W
23. ]
24. slope, intercept = calibrate_photodiode(cal_points)
25. # slope ≈ 5.0 (W/V), intercept ≈ 0.0
26. config.photodiode_slope = slope
27. config.photodiode_intercept = intercept

```

Step 4: Verification

```

# Verify calibration accuracy
def verify_calibration(slope, intercept, test_points):
    """Test calibration at intermediate points."""
    max_error = 0.0

    for voltage, expected_power in test_points:
        calculated_power = slope * voltage + intercept
        error = abs(calculated_power - expected_power) / expected_power

        if error > max_error:
            max_error = error

    return max_error

# Test at 3.75W (mid-point)
test_points = [(0.75, 3.75)]
max_error = verify_calibration(slope, intercept, test_points)

assert max_error < 0.05 # Must be within ±5%

```

Step 5: Documentation

```

# Log calibration to database
calibration_record = {
    "component": "photodiode",
    "date": datetime.now().isoformat(),
    "technician_id": current_user.id,
    "reference_equipment": "Ophir Nova II P/N 12345 (Cert: NIST-2025-001)",
    "calibration_points": cal_points,
    "slope": slope,
    "intercept": intercept,
    "max_error_percent": max_error * 100,
    "pass_fail": "PASS" if max_error < 0.05 else "FAIL",
    "notes": "Annual photodiode calibration"
}

db.save_calibration_record(calibration_record)

```

Acceptance Criteria: - [DONE] All 5 points within ±5% of certified meter - [DONE] Linear fit $R^2 > 0.99$ - [DONE] Verification point within ±5%

Actuator Calibration

Purpose

Verify actuator position accuracy (±10 μm tolerance).

Safety Critical: NO - Position errors do not directly harm patient

Equipment Required

1. **Digital micrometer** (0-25mm, $\pm 1\text{ }\mu\text{m}$ accuracy, NIST traceable)
2. **Gauge block** (10.000 mm $\pm 0.5\text{ }\mu\text{m}$, calibrated)
3. **TOSCA actuator** (Xeryon linear stage)

Procedure

Step 1: Homing

```
# Home actuator to establish reference
actuator.home()
assert actuator.is_homed
```

Step 2: Position Verification

Test at 5 positions across range:

Commanded Position (mm)	Micrometer Reading (mm)	Error (μm)	Acceptable?
0.000	0.000 ± 0.001	0	[DONE] $<10\text{ }\mu\text{m}$
5.000	5.002 ± 0.001	+2	[DONE] $<10\text{ }\mu\text{m}$
10.000	10.001 ± 0.001	+1	[DONE] $<10\text{ }\mu\text{m}$
15.000	14.998 ± 0.001	-2	[DONE] $<10\text{ }\mu\text{m}$
20.000	20.003 ± 0.001	+3	[DONE] $<10\text{ }\mu\text{m}$

Step 3: Calibration Adjustment

```
# If systematic error detected, apply offset
def adjust_actuator_offset(measured_errors: List[float]) -> float:
    """Calculate position offset from calibration errors."""
    mean_error = np.mean(measured_errors) # Average error in  $\mu\text{m}$ 

    if abs(mean_error) > 5.0: # If  $>5\text{ }\mu\text{m}$  average error
        return -mean_error # Correction offset
    else:
        return 0.0 # No correction needed

# Example
errors = [0, +2, +1, -2, +3] #  $\mu\text{m}$ 
offset = adjust_actuator_offset(errors) # offset =  $-0.8\text{ }\mu\text{m}$ 

config.actuator_position_offset_um = offset
```

Acceptance Criteria: - [DONE] All positions within $\pm 10\text{ }\mu\text{m}$ - [DONE] No systematic drift (error should average ~ 0)

Camera Calibration

Purpose

Verify pixel-to-micron calibration for ring detection.

Safety Critical: NO - Measurement only, not treatment control

Equipment Required

1. **Calibration target** (1951 USAF resolution chart or grid)
2. **Known dimensions** (e.g., 1.00 mm grid spacing)
3. **TOSCA camera** (Allied Vision)

Procedure

Step 1: Image Calibration Target

```
# Capture image of calibration target
frame = camera.capture_frame()
```

Step 2: Measure Pixel Spacing

```
# Detect grid lines in image
grid_spacing_pixels = detect_grid_spacing(frame) # e.g., 100 pixels

# Known physical spacing
grid_spacing_mm = 1.00 # 1 mm grid

# Calculate microns per pixel
um_per_pixel = (grid_spacing_mm * 1000) / grid_spacing_pixels
# Example: 1000 µm / 100 pixels = 10 µm/pixel

config.camera_um_per_pixel = um_per_pixel
```

Acceptance Criteria: - [DONE] Calculated µm/pixel within ±1% of expected - [DONE] Distortion <2% across image

Power Measurement Verification

Purpose

Verify laser power measurement accuracy during treatment.

Safety Critical: YES - Power limits enforced by this measurement

Procedure (Weekly)

```
# Quick verification using internal photodiode
def weekly_power_check():
    """Weekly verification of power measurement."""

    # Set laser to 5.0W
    laser.set_power(5.0)
    time.sleep(1.0) # Allow stabilization

    # Read photodiode
    measured_power = photodiode.get_power()

    # Verify within ±10%
    error = abs(measured_power - 5.0) / 5.0

    if error > 0.10:
        logger.error(f"Power measurement out of spec: {error*100:.1f}%")
        return False

    logger.info(f"Power check passed: {measured_power:.2f}W")
    return True
```

Acceptance Criteria: - [DONE] Measured power within ±10% of commanded

Calibration Schedule

Routine Calibration

Component	Frequency	Duration	Technician Required
Photodiode	Monthly	1 hour	Qualified technician

Actuator	Quarterly	30 min	Qualified technician
Camera	Annually	1 hour	Qualified technician
Power Check	Weekly	5 min	Operator
Watchdog	Annually	30 min	Qualified technician

Calibration Records

Required Documentation: 1. Calibration date and time 2. Technician name and ID 3. Reference equipment used (with certification numbers) 4. Calibration data (measurements, calculations) 5. Pass/Fail status 6. Next calibration due date

Database Schema:

```
CREATE TABLE calibration_records (  
  id INTEGER PRIMARY KEY AUTOINCREMENT,  
  component TEXT NOT NULL,  
  calibration_date TEXT NOT NULL,  
  technician_id INTEGER NOT NULL,  
  reference_equipment TEXT NOT NULL,  
  calibration_data TEXT, -- JSON  
  pass_fail TEXT NOT NULL,  
  next_due_date TEXT NOT NULL,  
  notes TEXT,  
  FOREIGN KEY (technician_id) REFERENCES technicians(id)  
);
```

Documentation & Audit Trail

FDA Requirements

21 CFR Part 820.72 - Inspection, Measuring, and Test Equipment:

Each manufacturer shall ensure that all inspection, measuring, and test equipment... is routinely calibrated, inspected, checked, and maintained.

TOSCA Compliance: - [DONE] Written calibration procedures (this document) - [DONE] Calibration schedule defined - [DONE] NIST traceable standards - [DONE] Documented calibration records - [DONE] Out-of-calibration handling procedures

Calibration Certificate Example

TOSCA CALIBRATION CERTIFICATE

Component: Photodiode Power Measurement

Serial Number: TOSCA-PD-001

Calibration Date: 2025-10-26

Technician: John Smith (Tech ID: 005)

Reference Equipment:

- Ophir Nova II Power Meter (S/N 12345)
- NIST Calibration Certificate: NIST-2025-001
- Expiration: 2026-10-15

Calibration Results:

- Test Points: 5 (0W, 2.5W, 5W, 7.5W, 10W)
- Maximum Error: 3.2% at 7.5W
- Calibration Curve: $P = 5.02V + 0.01$
- R^2 : 0.9997

Status: PASS (all points within ±5%)

Next Calibration Due: 2025-11-26

Signature: _____ Date: _____

References

Standards

- **ISO 10012:** Measurement management systems
- **ISO 13485:** Medical devices — Quality management
- **FDA 21 CFR Part 820:** Quality System Regulation
- **NIST Handbook 150:** NVLAP Procedures and General Requirements

Document Owner: Quality Assurance Manager **Last Updated:** 2025-10-26 **Next Review:** Before Phase 6 (calibration implementation) **Status:** Planning - Procedures defined, implementation pending