


# TOSCA Protocol Builder - Enhanced Action-Based Design

- [Architecture Diagrams](#)
  - [Figure 1: TOSCA Treatment Workflow Sequence](#)
- [Overview](#)
- [Enhanced Protocol Data Model](#)
  - [Action Types](#)
  - [JSON Schema](#)
  - [Action Parameter Schemas](#)
- [Example Protocols](#)
  - [Simple Treatment](#)
  - [Complex with Loops](#)
  - [Ramp with Actuator Movement](#)
- [GUI Design](#)
  - [Protocol Builder Tab Layout](#)
  - [Action Editor Dialog](#)
- [Protocol Execution](#)
  - [Execution Engine \(src/core/protocol\\_engine.py\)](#)
  - [Execution Recording](#)
- [Safety Checks](#)
- [File Storage](#)
- [Implementation Phases](#)

## Architecture Diagrams

**Figure 1: TOSCA Treatment Workflow Sequence**



TOSCA Treatment Workflow Sequence

Last Updated: 2025-11-04

Document Version: 2.0 Date: 2025-10-22 Supersedes: Section of 04\_treatment\_protocols.md

## Overview

The Protocol Builder provides a visual interface for creating action-based treatment protocols. Unlike the step-based model, this uses granular actions that can be sequenced, looped, and recorded.

## Enhanced Protocol Data Model

### Action Types

1. **SetLaserPower** - Instantly set laser to specific power
2. **RampLaserPower** - Ramp power from A to B over time
3. **MoveActuator** - Move to position at specified speed
4. **Wait** - Delay for specified duration
5. **Loop** - Repeat a block of actions N times or infinitely

### JSON Schema

```

{
  "$schema": "http://json-schema.org/draft-07/schema#",
  "type": "object",
  "required": ["protocol_name", "version", "actions"],
  "properties": {
  
```

```

"protocol_name": {"type": "string"},
"version": {"type": "string"},
"description": {"type": "string"},
"created_date": {"type": "string", "format": "date-time"},
"author": {"type": "string"},
"safety_limits": {
  "type": "object",
  "properties": {
    "max_power_watts": {"type": "number"},
    "max_duration_seconds": {"type": "number"},
    "min_actuator_position_um": {"type": "number"},
    "max_actuator_position_um": {"type": "number"}
  }
},
"actions": {
  "type": "array",
  "items": {
    "type": "object",
    "required": ["action_id", "action_type"],
    "properties": {
      "action_id": {"type": "integer"},
      "action_type": {
        "type": "string",
        "enum": ["SetLaserPower", "RampLaserPower", "MoveActuator", "Wait", "Loop"]
      },
      "parameters": {"type": "object"},
      "notes": {"type": "string"}
    }
  }
}
}

```

## Action Parameter Schemas

### SetLaserPower

```

{
  "action_type": "SetLaserPower",
  "parameters": {
    "power_watts": 5.0
  }
}

```

### RampLaserPower

```

{
  "action_type": "RampLaserPower",
  "parameters": {
    "start_power_watts": 2.0,
    "end_power_watts": 8.0,
    "duration_seconds": 30.0,
    "ramp_type": "linear" // "linear", "logarithmic", "exponential"
  }
}

```

### MoveActuator

```

{
  "action_type": "MoveActuator",
  "parameters": {
    "target_position_um": 2500,
    "speed_um_per_sec": 100
  }
}

```

### Wait

```

{
  "action_type": "Wait",
  "parameters": {
    "duration_seconds": 5.0
  }
}

```

### Loop

```

{
  "action_type": "Loop",
  "parameters": {
    "repeat_count": 3, // or -1 for infinite
    "actions": [
      // Nested actions go here
    ]
  }
}

```

## Example Protocols

### Simple Treatment

```

{
  "protocol_name": "Simple_5W_60s",
  "version": "1.0.0",
  "description": "Set 5W for 60 seconds",
  "created_date": "2025-10-22T00:00:00Z",
  "author": "User",
  "safety_limits": {
    "max_power_watts": 10.0,

```

```

    "max_duration_seconds": 300
  },
  "actions": [
    {
      "action_id": 1,
      "action_type": "SetLaserPower",
      "parameters": {"power_watts": 5.0},
      "notes": "Set to treatment power"
    },
    {
      "action_id": 2,
      "action_type": "Wait",
      "parameters": {"duration_seconds": 60.0},
      "notes": "Treatment duration"
    },
    {
      "action_id": 3,
      "action_type": "SetLaserPower",
      "parameters": {"power_watts": 0.0},
      "notes": "Turn off laser"
    }
  ]
}

```

## Complex with Loops

```

{
  "protocol_name": "Pulsed_Treatment_3x",
  "version": "1.0.0",
  "description": "3 pulses of 5W, 10s each, with 5s gaps",
  "created_date": "2025-10-22T00:00:00Z",
  "actions": [
    {
      "action_id": 1,
      "action_type": "MoveActuator",
      "parameters": {
        "target_position_um": 2000,
        "speed_um_per_sec": 50
      }
    },
    {
      "action_id": 2,
      "action_type": "Loop",
      "parameters": {
        "repeat_count": 3,
        "actions": [
          {
            "action_id": 3,
            "action_type": "SetLaserPower",
            "parameters": {"power_watts": 5.0}
          },
          {
            "action_id": 4,
            "action_type": "Wait",
            "parameters": {"duration_seconds": 10.0}
          },
          {
            "action_id": 5,
            "action_type": "SetLaserPower",
            "parameters": {"power_watts": 0.0}
          },
          {
            "action_id": 6,
            "action_type": "Wait",
            "parameters": {"duration_seconds": 5.0}
          }
        ]
      }
    }
  ]
}

```

## Ramp with Actuator Movement

```

{
  "protocol_name": "Ramp_With_Ring_Adjust",
  "version": "1.0.0",
  "description": "Ramp power while adjusting ring size",
  "created_date": "2025-10-22T00:00:00Z",
  "actions": [
    {
      "action_id": 1,
      "action_type": "MoveActuator",
      "parameters": {
        "target_position_um": 3000,
        "speed_um_per_sec": 100
      },
      "notes": "Start at 3mm ring"
    },
    {
      "action_id": 2,
      "action_type": "RampLaserPower",
      "parameters": {
        "start_power_watts": 2.0,
        "end_power_watts": 8.0,
        "duration_seconds": 60.0,
        "ramp_type": "linear"
      },
      "notes": "Ramp power over 60 seconds"
    },
    {
      "action_id": 3,
      "action_type": "MoveActuator",
      "parameters": {
        "target_position_um": 2000,

```

```
        "speed_um_per_sec": 50
    },
    "notes": "Reduce to 2mm ring"
},
{
    "action_id": 4,
    "action_type": "Wait",
    "parameters": {"duration_seconds": 30.0}
},
{
    "action_id": 5,
    "action_type": "SetLaserPower",
    "parameters": {"power_watts": 0.0}
}
]
}
```

## GUI Design

### Protocol Builder Tab Layout

```
1. **Protocol Builder [?] [ø]**
2. **Protocol** - [Unnamed Protocol] [Open] [Save] [Save As]
3. **Description** - [ ]
4. **Action Sequence**
5. **# Type Parameters [][ ][x]**
6. **1 SetLaserPower 5.0 W [][ ][x]**
7. **2 Wait 10.0 s [][ ][x]**
8. **3 MoveActuator 2500 µm @ 100 µ/s [][ ][x]**
9. **4 Loop (3x) ▼ 4 actions [][ ][x]**
10. **RampLaserPower 28W, 30s**
11. **Wait 5s**
12. **SetLaserPower 0W**
13. **Wait 2s**
14. **[+ Add Action ▼]**
15. *** Set Laser Power**
16. *** Ramp Laser Power**
17. *** Move Actuator**
18. *** Wait**
19. *** Loop**
20. **Total Duration** - ~120s Max Power: 8.0W
21. **[Execute]**
22. **Status** - Ready Safety: [DONE] OK [Execute & Rec]
```

### Action Editor Dialog

When user clicks [🔍] or [+ Add Action], show dialog:

```
Edit Action: RampLaserPower

Start Power: [2.0 ] W
End Power: [8.0 ] W
Duration: [30.0 ] seconds
Ramp Type: [Linear ▼]
            • Linear
            • Logarithmic
            • Exponential

Notes: [Warm-up phase_____]

[Preview Curve]

[Cancel] [OK]
```

## Protocol Execution

### Execution Engine (src/core/protocol\_engine.py)

```
class ProtocolEngine:
    """Execute action-based protocols with real-time control"""

    def __init__(self, laser_controller, actuator_controller):
        self.laser = laser_controller
        self.actuator = actuator_controller
        self.is_running = False
        self.current_action_id = None

    async def execute_protocol(self, protocol: dict, record: bool = False):
        """Execute protocol actions sequentially"""

    async def execute_action(self, action: dict):
        """Execute single action"""

    async def execute_loop(self, loop_action: dict):
        """Execute looped actions"""

    def pause(self):
        """Pause execution"""

    def resume(self):
        """Resume execution"""

    def stop(self):
        """Stop and reset"""
```

### Execution Recording

Record to database table `protocol_executions`: - `execution_id` - `protocol_name` - `start_time` - `end_time` - `subject_id` (optional) - `session_id` - `execution_log` (JSON: action timestamps, actual values) - `video_file_path`

## Safety Checks

Before execution: 1. Validate all parameters against `safety_limits` 2. Check footpedal status 3. Verify subject selected 4. Confirm actuator homed and in range 5. Verify camera feed active

During execution: - Monitor safety interlocks continuously - Log every action with timestamp - Allow emergency stop at any time - Validate each action before execution

## File Storage

Protocols saved to: `data/protocols/` - Format: JSON - Naming: `{protocol_name}_v{version}.json` - Include metadata: creation date, author, modifications

## Implementation Phases

1. **Phase 1: Data Model & Core Engine**
  - Define action classes
  - Implement `protocol_engine.py`
  - Add JSON validation
2. **Phase 2: Basic GUI**
  - Create `protocol_builder_widget.py`
  - Action list table
  - Simple add/remove/edit
3. **Phase 3: Advanced Features**
  - Loop support with nesting
  - Drag-and-drop reordering
  - Protocol templates library
4. **Phase 4: Execution & Recording**
  - Real-time execution
  - Progress visualization
  - Execution recording to database
5. **Phase 5: Polish**
  - Parameter validation UI
  - Curve preview for ramps
  - Import/export protocols

---

**Next Steps:** 1. Review and approve this design 2. Implement protocol data classes in `src/core/protocol.py` 3. Create `protocol_builder_widget.py` UI skeleton 4. Begin Phase 1 implementation