

TOSCA Laser Control System - Image Processing Pipeline

- [Architecture Diagrams](#)
- [Overview](#)
- [Camera System](#)
 - [Hardware](#)
 - [Camera Controller](#)
- [Ring Detection](#)
 - [Algorithm: Hough Circle Transform](#)
- [Focus Measurement](#)
 - [Algorithm: Laplacian Variance](#)
- [Video Recording](#)
 - [VideoRecorder Class](#)
- [Image Processing Pipeline Integration](#)
 - [Complete Frame Processor](#)
- [UI Integration](#)
 - [Video Display Widget \(PyQt6\)](#)
- [Calibration Procedures](#)
 - [1. Ring Size Calibration](#)
 - [2. Focus Threshold Calibration](#)
- [Performance Optimization](#)
 - [Tips for Real-Time Processing](#)

Architecture Diagrams

Document Version: 1.0 **Date:** 2025-10-15

Overview

The image processing pipeline handles camera feed processing for: - **Ring detection** - Locating laser ring position and size - **Focus measurement** - Quantifying image sharpness - **Video recording** - Capturing treatment sessions - **Live display** - Real-time feed with overlays

All processing is done using OpenCV (cv2) and NumPy.

Camera System

Hardware

- **SDK:** VmbPy (Allied Vision Vimba Python SDK)
- **Interface:** USB 3.0 / GigE
- **Resolution:** 1920x1080 (configurable)
- **Frame rate:** 30 FPS (configurable)
- **Color:** Monochrome or Color (prefer monochrome for laser applications)

Camera Controller

```
from vmbpy import *
import numpy as np
import threading
import queue

class CameraController:
```

```

"""
Wrapper for VmbPy camera control
"""

def __init__(self, camera_id=None):
    self.vmb = VmbSystem.get_instance()
    self.camera = None
    self.is_streaming = False
    self.frame_queue = queue.Queue(maxsize=10)
    self.frame_callback = None

    # Camera settings
    self.exposure_us = 10000 # 10ms default
    self.gain_db = 0
    self.resolution = (1920, 1080)
    self.fps = 30

def connect(self, camera_id=None) -> tuple[bool, str]:
    """Connect to camera"""
    try:
        with self.vmb:
            cameras = self.vmb.get_all_cameras()
            if not cameras:
                return (False, "No cameras found")

            # Use specified camera or first available
            if camera_id:
                self.camera = self.vmb.get_camera_by_id(camera_id)
            else:
                self.camera = cameras[0]

            self.camera.open()

            # Configure camera
            self._configure_camera()

            return (True, f"Connected to {self.camera.get_id()}")

    except Exception as e:
        return (False, f"Camera connection error: {str(e)}")

def _configure_camera(self):
    """Set camera parameters"""
    # Set exposure
    self.camera.ExposureTime.set(self.exposure_us)

    # Set gain
    self.camera.Gain.set(self.gain_db)

    # Set resolution (if camera supports)
    try:
        self.camera.Width.set(self.resolution[0])
        self.camera.Height.set(self.resolution[1])
    except:
        pass # Use camera default

    # Set frame rate
    try:
        self.camera.AcquisitionFrameRate.set(self.fps)
    except:
        pass

def start_streaming(self):
    """Start continuous frame acquisition"""
    if not self.camera:
        return

    self.is_streaming = True
    self.camera.start_streaming(self._frame_handler)

def stop_streaming(self):
    """Stop frame acquisition"""
    if self.camera and self.is_streaming:
        self.camera.stop_streaming()

```

```

        self.is_streaming = False

def _frame_handler(self, cam, stream, frame):
    """
    Callback for each frame from camera

    Args:
        frame: VmbFrame object
    """
    try:
        # Convert to numpy array
        img_array = frame.as_numpy_ndarray()

        # Add to queue (non-blocking)
        if not self.frame_queue.full():
            self.frame_queue.put(img_array, block=False)

        # Call external callback if set
        if self.frame_callback:
            self.frame_callback(img_array)

    except queue.Full:
        pass # Drop frame if queue full
    except Exception as e:
        logger.error(f"Frame handler error: {e}")

def get_latest_frame(self) -> np.ndarray:
    """Get most recent frame from queue"""
    try:
        # Clear queue and get latest
        frame = None
        while not self.frame_queue.empty():
            frame = self.frame_queue.get_nowait()
        return frame
    except queue.Empty:
        return None

def disconnect(self):
    """Disconnect camera"""
    self.stop_streaming()
    if self.camera:
        self.camera.close()
        self.camera = None

```

Ring Detection

Algorithm: Hough Circle Transform

The laser produces a circular ring that must be detected for alignment verification.

```

import cv2
import numpy as np
from dataclasses import dataclass

@dataclass
class DetectedRing:
    """Result of ring detection"""
    center_x: int
    center_y: int
    radius: int
    confidence: float # 0.0 to 1.0
    detected: bool

class RingDetector:
    """
    Detect circular laser ring in camera images
    """

    def __init__(self):
        # Detection parameters (tunable)
        self.dp = 1 # Inverse ratio of accumulator resolution
        self.min_dist = 100 # Minimum distance between circle centers

```

```

self.param1 = 50 # Canny edge detector high threshold
self.param2 = 30 # Accumulator threshold for circle centers
self.min_radius = 30 # Minimum circle radius (pixels)
self.max_radius = 400 # Maximum circle radius (pixels)

# Preprocessing parameters
self.blur_kernel = (9, 9)
self.blur_sigma = 2

# Region of interest (optional - speeds up detection)
self.roi = None # (x, y, width, height) or None for full frame

def detect(self, frame: np.ndarray) -> DetectedRing:
    """
    Detect laser ring in frame

    Args:
        frame: Input image (color or grayscale)

    Returns:
        DetectedRing object
    """
    # Convert to grayscale if needed
    if len(frame.shape) == 3:
        gray = cv2.cvtColor(frame, cv2.COLOR_BGR2GRAY)
    else:
        gray = frame.copy()

    # Apply ROI if set
    if self.roi:
        x, y, w, h = self.roi
        gray = gray[y:y+h, x:x+w]
    else:
        x, y = 0, 0

    # Preprocessing
    preprocessed = self._preprocess(gray)

    # Detect circles using Hough Transform
    circles = cv2.HoughCircles(
        preprocessed,
        cv2.HOUGH_GRADIENT,
        dp=self.dp,
        minDist=self.min_dist,
        param1=self.param1,
        param2=self.param2,
        minRadius=self.min_radius,
        maxRadius=self.max_radius
    )

    if circles is None or len(circles) == 0:
        return DetectedRing(0, 0, 0, 0.0, False)

    # Select best circle
    best_circle = self._select_best_circle(circles[0], preprocessed)

    # Adjust coordinates for ROI offset
    return DetectedRing(
        center_x=int(best_circle[0]) + x,
        center_y=int(best_circle[1]) + y,
        radius=int(best_circle[2]),
        confidence=self._calculate_confidence(best_circle, preprocessed),
        detected=True
    )

def _preprocess(self, gray: np.ndarray) -> np.ndarray:
    """Preprocess image for better circle detection"""

    # 1. Gaussian blur to reduce noise
    blurred = cv2.GaussianBlur(gray, self.blur_kernel, self.blur_sigma)

    # 2. Enhance contrast (optional - adjust if needed)
    # clahe = cv2.createCLAHE(clipLimit=2.0, tileGridSize=(8,8))
    # enhanced = clahe.apply(blurred)

```

```

    return blurred

def _select_best_circle(self, circles: np.ndarray, image: np.ndarray):
    """
    Select best circle from multiple detections

    Criteria:
    - Brightness of ring perimeter
    - Darkness of interior
    - Circularity
    """
    if len(circles) == 1:
        return circles[0]

    best_score = -1
    best_circle = circles[0]

    for circle in circles:
        score = self._score_circle(circle, image)
        if score > best_score:
            best_score = score
            best_circle = circle

    return best_circle

def _score_circle(self, circle, image: np.ndarray) -> float:
    """
    Score circle quality

    Higher score = better match to laser ring
    """
    cx, cy, r = circle

    # Create mask for ring perimeter
    mask_outer = np.zeros(image.shape, dtype=np.uint8)
    cv2.circle(mask_outer, (int(cx), int(cy)), int(r), 255, thickness=3)

    # Create mask for interior
    mask_inner = np.zeros(image.shape, dtype=np.uint8)
    cv2.circle(mask_inner, (int(cx), int(cy)), int(r - 10), 255, thickness=-1)

    # Calculate mean intensity of ring vs interior
    ring_intensity = cv2.mean(image, mask=mask_outer)[0]
    inner_intensity = cv2.mean(image, mask=mask_inner)[0]

    # Good laser ring: bright perimeter, darker interior
    contrast = ring_intensity - inner_intensity

    # Normalize score
    score = contrast / 255.0

    return max(0.0, min(1.0, score))

def _calculate_confidence(self, circle, image: np.ndarray) -> float:
    """Calculate detection confidence (0.0 to 1.0)"""
    return self._score_circle(circle, image)

def draw_overlay(self, frame: np.ndarray, ring: DetectedRing) -> np.ndarray:
    """
    Draw detected ring on frame

    Returns:
    Frame with overlay
    """
    overlay = frame.copy()

    if not ring.detected:
        # No ring detected - draw warning
        cv2.putText(overlay, "RING NOT DETECTED", (50, 50),
                    cv2.FONT_HERSHEY_SIMPLEX, 1, (0, 0, 255), 2)
        return overlay

    # Draw circle

```

```

color = (0, 255, 0) if ring.confidence > 0.5 else (0, 165, 255) # Green if confident, orange if
not
cv2.circle(overlay, (ring.center_x, ring.center_y), ring.radius, color, 2)

# Draw center crosshair
cv2.drawMarker(overlay, (ring.center_x, ring.center_y), color,
               cv2.MARKER_CROSS, 20, 2)

# Draw info text
info_text = f"Ring: {ring.radius}px ({ring.confidence:.2f})"
cv2.putText(overlay, info_text, (ring.center_x - 50, ring.center_y - ring.radius - 10),
            cv2.FONT_HERSHEY_SIMPLEX, 0.6, color, 2)

return overlay

def calibrate_ring_size(self, ring: DetectedRing, actuator_position_um: float) -> dict:
    """
    Create calibration point mapping ring pixel radius to physical size

    Returns:
        Calibration data point
    """
    return {
        'actuator_position_um': actuator_position_um,
        'ring_radius_pixels': ring.radius,
        'center_x': ring.center_x,
        'center_y': ring.center_y,
        'confidence': ring.confidence
    }

```

Focus Measurement

Algorithm: Laplacian Variance

Focus quality is measured using the variance of the Laplacian operator - sharper images have higher variance.

```

class FocusAnalyzer:
    """
    Measure image focus quality
    """

    def __init__(self):
        # Calibration - these values depend on your camera/optics
        self.min_score = 0 # Completely out of focus
        self.max_score = 1000 # Perfectly in focus
        self.good_focus_threshold = 300 # Minimum for "acceptable" focus

    def calculate_focus_score(self, frame: np.ndarray) -> float:
        """
        Calculate focus quality score

        Args:
            frame: Input image (color or grayscale)

        Returns:
            Focus score (higher = better focus)
        """
        # Convert to grayscale if needed
        if len(frame.shape) == 3:
            gray = cv2.cvtColor(frame, cv2.COLOR_BGR2GRAY)
        else:
            gray = frame

        # Calculate Laplacian
        laplacian = cv2.Laplacian(gray, cv2.CV_64F)

        # Calculate variance
        variance = laplacian.var()

        return variance

```

```

def calculate_focus_percent(self, frame: np.ndarray) -> float:
    """
    Calculate focus as percentage (0-100)

    Returns:
        Focus percentage (100 = excellent focus, 0 = very poor)
    """
    score = self.calculate_focus_score(frame)

    # Normalize to 0-100 scale
    percent = ((score - self.min_score) / (self.max_score - self.min_score)) * 100
    percent = max(0, min(100, percent))

    return percent

def is_in_focus(self, frame: np.ndarray) -> tuple[bool, float]:
    """
    Determine if image is acceptably in focus

    Returns:
        (is_focused, focus_score)
    """
    score = self.calculate_focus_score(frame)
    is_focused = score >= self.good_focus_threshold

    return (is_focused, score)

def draw_focus_indicator(self, frame: np.ndarray) -> np.ndarray:
    """
    Draw focus quality indicator on frame

    Returns:
        Frame with focus indicator overlay
    """
    overlay = frame.copy()

    # Calculate focus
    is_focused, score = self.is_in_focus(frame)
    focus_percent = self.calculate_focus_percent(frame)

    # Draw focus bar (top-left corner)
    bar_x, bar_y = 20, 20
    bar_width, bar_height = 200, 30

    # Background
    cv2.rectangle(overlay, (bar_x, bar_y), (bar_x + bar_width, bar_y + bar_height),
                  (100, 100, 100), -1)

    # Fill based on focus percentage
    fill_width = int(bar_width * (focus_percent / 100))
    color = (0, 255, 0) if is_focused else (0, 165, 255) # Green if focused, orange if not
    cv2.rectangle(overlay, (bar_x, bar_y), (bar_x + fill_width, bar_y + bar_height),
                  color, -1)

    # Border
    cv2.rectangle(overlay, (bar_x, bar_y), (bar_x + bar_width, bar_y + bar_height),
                  (255, 255, 255), 2)

    # Text
    text = f"Focus: {focus_percent:.0f}%"
    cv2.putText(overlay, text, (bar_x + 5, bar_y + 20),
                cv2.FONT_HERSHEY_SIMPLEX, 0.6, (255, 255, 255), 2)

    return overlay

def calibrate_focus_thresholds(self, focused_frames: list, unfocused_frames: list):
    """
    Calibrate focus thresholds based on sample images

    Args:
        focused_frames: List of frames that are in good focus
        unfocused_frames: List of frames that are out of focus
    """
    # Calculate scores for focused frames

```

```

focused_scores = [self.calculate_focus_score(frame) for frame in focused_frames]
min_focused = min(focused_scores)
max_focused = max(focused_scores)

# Calculate scores for unfocused frames
unfocused_scores = [self.calculate_focus_score(frame) for frame in unfocused_frames]
max_unfocused = max(unfocused_scores)

# Set threshold between max unfocused and min focused
self.good_focus_threshold = (max_unfocused + min_focused) / 2
self.max_score = max_focused * 1.5 # Leave headroom

print(f"Calibration results:")
print(f"  Focused range: {min_focused:.1f} - {max_focused:.1f}")
print(f"  Unfocused range: {min(unfocused_scores):.1f} - {max_unfocused:.1f}")
print(f"  Threshold set to: {self.good_focus_threshold:.1f}")

```

Video Recording

VideoRecorder Class

```

class VideoRecorder:
    """
    Record video from camera feed
    """

    def __init__(self, session_folder: str):
        self.session_folder = session_folder
        self.video_writer = None
        self.is_recording = False
        self.frame_count = 0
        self.start_time = None

        # Video settings
        self.fourcc = cv2.VideoWriter_fourcc(*'XVID') # or 'mp4v', 'H264'
        self.fps = 30
        self.resolution = (1920, 1080)

    def start_recording(self, filename: str = None) -> tuple[bool, str]:
        """
        Start video recording

        Args:
            filename: Optional custom filename (default: video.avi)

        Returns:
            (success, filepath)
        """
        if self.is_recording:
            return (False, "Already recording")

        # Generate filename
        if not filename:
            timestamp = datetime.now().strftime("%Y%m%d_%H%M%S")
            filename = f"treatment_{timestamp}.avi"

        filepath = os.path.join(self.session_folder, filename)

        # Create VideoWriter
        try:
            self.video_writer = cv2.VideoWriter(
                filepath,
                self.fourcc,
                self.fps,
                self.resolution
            )

            if not self.video_writer.isOpened():
                return (False, "Failed to open video writer")

            self.is_recording = True
            self.frame_count = 0

```



```

        self.start_time = time.time()

        return (True, filepath)

    except Exception as e:
        return (False, f"Recording start error: {str(e)}")

def write_frame(self, frame: np.ndarray):
    """
    Write frame to video file

    Args:
        frame: Image frame (must match resolution)
    """
    if not self.is_recording or self.video_writer is None:
        return

    # Resize if needed
    if frame.shape[:2] != (self.resolution[1], self.resolution[0]):
        frame = cv2.resize(frame, self.resolution)

    # Write frame
    self.video_writer.write(frame)
    self.frame_count += 1

def stop_recording(self) -> dict:
    """
    Stop recording and finalize video

    Returns:
        Recording metadata dict
    """
    if not self.is_recording:
        return {}

    duration = time.time() - self.start_time

    # Release video writer
    if self.video_writer:
        self.video_writer.release()
        self.video_writer = None

    self.is_recording = False

    metadata = {
        'frame_count': self.frame_count,
        'duration_seconds': duration,
        'fps': self.fps,
        'resolution': f"{self.resolution[0]}x{self.resolution[1]}"
    }

    return metadata

def capture_snapshot(self, frame: np.ndarray, label: str = "") -> str:
    """
    Save single frame as image

    Args:
        frame: Image to save
        label: Optional label for filename

    Returns:
        Filepath of saved image
    """
    timestamp = datetime.now().strftime("%Y%m%d_%H%M%S%f")
    if label:
        filename = f"snapshot_{label}_{timestamp}.png"
    else:
        filename = f"snapshot_{timestamp}.png"

    filepath = os.path.join(self.session_folder, "snapshots", filename)

    # Create directory if needed
    os.makedirs(os.path.dirname(filepath), exist_ok=True)

```

```

# Save image
cv2.imwrite(filepath, frame)

return filepath

```

Image Processing Pipeline Integration

Complete Frame Processor

```

class FrameProcessor:
    """
    Complete image processing pipeline
    """

    def __init__(self, session_folder: str):
        self.ring_detector = RingDetector()
        self.focus_analyzer = FocusAnalyzer()
        self.video_recorder = VideoRecorder(session_folder)

        self.enable_ring_detection = True
        self.enable_focus_display = True
        self.enable_recording = False

    def process_frame(self, frame: np.ndarray) -> tuple[np.ndarray, dict]:
        """
        Process camera frame

        Args:
            frame: Raw camera frame

        Returns:
            (processed_frame, metadata)
        """
        processed = frame.copy()
        metadata = {}

        # Ring detection
        if self.enable_ring_detection:
            ring = self.ring_detector.detect(frame)
            processed = self.ring_detector.draw_overlay(processed, ring)
            metadata['ring'] = {
                'detected': ring.detected,
                'center_x': ring.center_x,
                'center_y': ring.center_y,
                'radius': ring.radius,
                'confidence': ring.confidence
            }

        # Focus measurement
        if self.enable_focus_display:
            is_focused, focus_score = self.focus_analyzer.is_in_focus(frame)
            processed = self.focus_analyzer.draw_focus_indicator(processed)
            metadata['focus'] = {
                'is_focused': is_focused,
                'score': focus_score,
                'percent': self.focus_analyzer.calculate_focus_percent(frame)
            }

        # Record frame if recording active
        if self.enable_recording and self.video_recorder.is_recording:
            self.video_recorder.write_frame(processed)

        return (processed, metadata)

    def start_recording(self):
        """Start video recording"""
        self.enable_recording = True
        return self.video_recorder.start_recording()

    def stop_recording(self):
        """Stop video recording"""

```

```

self.enable_recording = False
return self.video_recorder.stop_recording()

```

UI Integration

Video Display Widget (PyQt6)

```

from PyQt6.QtWidgets import QLabel, QWidget, QVBoxLayout
from PyQt6.QtGui import QImage, QPixmap
from PyQt6.QtCore import QTimer, pyqtSignal
import cv2

class VideoDisplayWidget(QWidget):
    """
    PyQt6 widget for displaying live camera feed
    """

    frame_processed = pyqtSignal(dict) # Emit metadata

    def __init__(self, camera_controller, frame_processor):
        super().__init__()

        self.camera = camera_controller
        self.processor = frame_processor

        # UI setup
        self.image_label = QLabel()
        layout = QVBoxLayout()
        layout.addWidget(self.image_label)
        self.setLayout(layout)

        # Update timer
        self.timer = QTimer()
        self.timer.timeout.connect(self.update_frame)
        self.timer.start(33) # ~30 FPS

    def update_frame(self):
        """Update displayed frame"""
        # Get latest frame from camera
        frame = self.camera.get_latest_frame()
        if frame is None:
            return

        # Process frame
        processed, metadata = self.processor.process_frame(frame)

        # Emit metadata
        self.frame_processed.emit(metadata)

        # Convert to Qt format and display
        self.display_frame(processed)

    def display_frame(self, frame):
        """Convert OpenCV frame to Qt and display"""
        # Convert BGR to RGB
        rgb_frame = cv2.cvtColor(frame, cv2.COLOR_BGR2RGB)

        h, w, ch = rgb_frame.shape
        bytes_per_line = ch * w

        # Convert to QImage
        qt_image = QImage(rgb_frame.data, w, h, bytes_per_line, QImage.Format.Format_RGB888)

        # Display
        self.image_label.setPixmap(QPixmap.fromImage(qt_image))

```

Calibration Procedures

1. Ring Size Calibration

Goal: Map actuator position to physical ring diameter

```
def calibrate_ring_size():
    """Interactive ring size calibration"""

    calibration_points = []

    print("Ring Size Calibration")
    print("-----")

    # Move actuator through range
    positions_um = [0, 500, 1000, 1500, 2000, 2500]

    for position in positions_um:
        # Move actuator
        actuator.move_to_position(position)
        time.sleep(1) # Allow settling

        # Capture frame
        frame = camera.get_latest_frame()

        # Detect ring
        ring = ring_detector.detect(frame)

        if not ring.detected:
            print(f" Position {position}μm: RING NOT DETECTED")
            continue

        # Prompt user for actual physical measurement
        print(f" Position {position}μm: Detected {ring.radius}px")
        physical_size_mm = float(input(" Enter actual ring diameter (mm): "))

        calibration_points.append({
            'actuator_position_um': position,
            'ring_size_mm': physical_size_mm,
            'ring_radius_pixels': ring.radius
        })

    # Save calibration to database
    save_calibration('actuator_position_to_ring_size', calibration_points)

    print("Calibration complete!")
    return calibration_points
```

2. Focus Threshold Calibration

Goal: Determine focus score threshold for “in focus” determination

```
def calibrate_focus():
    """Interactive focus calibration"""

    print("Focus Calibration")
    print("-----")
    print("Adjust focus to BEST focus, then press Enter...")
    input()

    # Capture focused frames
    focused_frames = []
    for i in range(10):
        frame = camera.get_latest_frame()
        focused_frames.append(frame)
        time.sleep(0.1)

    print("Now adjust to OUT OF FOCUS, then press Enter...")
    input()

    # Capture unfocused frames
    unfocused_frames = []
    for i in range(10):
        frame = camera.get_latest_frame()
        unfocused_frames.append(frame)
        time.sleep(0.1)
```

```
# Calibrate thresholds
focus_analyzer.calibrate_focus_thresholds(focused_frames, unfocused_frames)

# Save calibration
calibration_data = {
    'good_focus_threshold': focus_analyzer.good_focus_threshold,
    'max_score': focus_analyzer.max_score
}
save_calibration('camera_focus_calibration', calibration_data)

print("Focus calibration complete!")
```

Performance Optimization

Tips for Real-Time Processing

1. **Use ROI:** Restrict ring detection to region of interest
 2. **Downscale:** Process smaller image for speed (e.g., 960x540 instead of 1920x1080)
 3. **Threading:** Process frames in separate thread from UI
 4. **Skip frames:** Don't process every frame if 30 FPS is too fast
 5. **GPU acceleration:** Use OpenCV with CUDA if available
-

Document Owner: Image Processing Engineer **Last Updated:** 2025-10-15