

TOSCA Implementation Status & Readiness Tiers

- [Architecture Diagrams](#)
- [Overview](#)
- [Tier Classification](#)
- [TIER 1: Lab/Experimentation \(\[DONE\] CURRENT STATUS\)](#)
 - [Implemented Features \(Phase 5\)](#)
- [TIER 2: Pre-Clinical Validation \(WARNING: PLANNED - Phase 6\)](#)
 - [Required Features \(Not Yet Implemented\)](#)
- [TIER 3: Clinical Trials \(\[FAILED\] NOT STARTED - Phase 7\)](#)
- [TIER 4: FDA Submission/Production \(\[FAILED\] NOT STARTED - Phase 8+\)](#)
- [Quick Reference: What Can I Do?](#)
 - [\[DONE\] Right Now \(Tier 1 - Current\)](#)
 - [WARNING: After Phase 6 \(Tier 2 - Pre-Clinical\)](#)
 - [After Phase 7 \(Tier 3 - Clinical Trials\)](#)
 - [After Phase 8+ \(Tier 4 - FDA/Production\)](#)
- [Development Roadmap](#)
- [Risk Assessment by Tier](#)
- [Summary](#)

Architecture Diagrams

Document Version: 1.2 **Last Updated:** 2025-11-05 **Current Version:** v0.9.13-alpha **Purpose:** Clear categorization of what's implemented vs what's required for different use cases

Overview

This document categorizes TOSCA features into **readiness tiers** based on what's needed for different stages of development and deployment.

Tier Classification

1. ****TIER 1**** - Lab/Experimentation (CURRENT - Phase 5)
 2. ****[DONE]** Implemented - Safe for R&D with test data**
 3. ****Need to add****
 4. ****TIER 2**** - Pre-Clinical Validation (Phase 6)
 5. ****WARNING**** - Partially Planned - Required for cadaver/bench testing
 6. ****Need to add****
 7. ****TIER 3**** - Clinical Trials (Phase 7)
 8. ****[FAILED]** Not Started - Required for human testing (IRB approval)**
 9. ****Need to add****
 10. ****TIER 4**** - FDA Submission/Production (Phase 8+)
 11. ****[FAILED]** Not Started - Required for commercial deployment**
-

TIER 1: Lab/Experimentation ([DONE] CURRENT STATUS)

What You Can Do: - [DONE] Internal R&D and proof-of-concept testing - [DONE] Algorithm development and optimization - [DONE] Bench testing with test patterns (no real patients) - [DONE] Software development and debugging - [DONE] Unit and integration testing - [DONE] Team training and familiarization

What You CANNOT Do: - [FAILED] Clinical trials with humans - [FAILED] Animal testing (requires validation) - [FAILED] Store real patient data (no encryption) - [FAILED] Production deployment - [FAILED] FDA submission

Implemented Features (Phase 5)

| Feature | Status | Document |
|--------------------------------------|----------------|---------------------------|
| Core Treatment Logic | | |
| Laser power control (5W) | [DONE] Working | 01_system_overview.md |
| Actuator position control (45mm) | [DONE] Working | 01_system_overview.md |
| Protocol execution engine | [DONE] Working | 06_protocol_builder.md |
| Treatment session management | [DONE] Working | - |
| Safety Systems | | |
| 7 independent interlocks | [DONE] Working | 03_safety_system.md |
| Selective shutdown (laser only) | [DONE] Working | SAFETY_SHUTDOWN_POLICY.md |
| Emergency stop handling | [DONE] Working | 03_safety_system.md |
| Safety watchdog (500ms heartbeat) | [DONE] Working | 07_safety_watchdog.md |
| GPIO monitoring (Arduino Nano) | [DONE] Working | 03_safety_system.md |
| Camera & Image Processing | | |
| Real-time streaming (30 FPS) | [DONE] Working | - |
| Allied Vision API integration | [DONE] Working | camera_controller.py |
| Thread-safe camera controls | [DONE] Working | camera_controller.py |
| Exposure/gain hardware feedback | [DONE] Working | camera_widget.py |
| Pixel format auto-detection | [DONE] Working | camera_controller.py |
| Ring detection | [DONE] Working | 05_image_processing.md |
| Focus measurement | [DONE] Working | 05_image_processing.md |
| Recording | | |
| Video recording (MP4, unencrypted) | [DONE] Working | 12_recording_manager.md |
| Data Management | | |
| SQLite database (unencrypted) | [DONE] Working | 02_database_schema.md |
| Session tracking | [DONE] Working | - |
| Treatment event logging | [DONE] Working | 11_event_logging.md |
| Two-tier logging (JSONL + DB) | [DONE] Working | 11_event_logging.md |
| Testing Infrastructure | | |
| MockHardwareBase pattern | [DONE] Working | 09_test_architecture.md |
| 30+ test files | [DONE] Working | 09_test_architecture.md |
| ~85% code coverage | [DONE] Working | - |
| Threading | | |
| Camera streaming thread | [DONE] Working | 10_concurrency_model.md |
| Safety watchdog timer | [DONE] Working | 10_concurrency_model.md |
| PyQt6 signal/slot communication | [DONE] Working | 10_concurrency_model.md |

Tier 1 Summary: 27 features implemented [DONE] (4 new camera features added Oct 2025)

TIER 2: Pre-Clinical Validation (WARNING: PLANNED - Phase 6)

What This Enables: - WARNING: Cadaver testing (dead tissue, no live patients) - WARNING: Bench validation with calibrated equipment - WARNING: Performance verification against specifications - WARNING: Safety system validation testing - WARNING: Pre-clinical animal testing (if applicable)

Still Cannot Do: - [FAILED] Human clinical trials (Tier 3 required) - [FAILED] FDA submission (Tier 4

required)

Required Features (Not Yet Implemented)

| Feature | Status | Priority | Document |
|---|------------------|----------|------------------------------|
| Security (Critical) | | | |
| Database encryption (SQLCipher AES-256) | WARNING: Planned | P0 | 08_security_architecture.md |
| Video file encryption (AES-256-GCM) | WARNING: Planned | P0 | 12_recording_manager.md |
| Configuration file encryption | WARNING: Planned | P1 | 08_security_architecture.md |
| Key derivation (PBKDF2) | WARNING: Planned | P0 | 08_security_architecture.md |
| Calibration (Critical) | | | |
| Photodiode calibration workflow | WARNING: Planned | P0 | 13_calibration_procedures.md |
| Actuator position calibration | WARNING: Planned | P0 | 13_calibration_procedures.md |
| Camera pixel calibration | WARNING: Planned | P1 | 13_calibration_procedures.md |
| Calibration database & records | WARNING: Planned | P0 | 13_calibration_procedures.md |
| NIST-traceable certificates | WARNING: Planned | P0 | 13_calibration_procedures.md |
| Audit Trail Integrity | | | |
| HMAC signatures per event | WARNING: Planned | P1 | 11_event_logging.md |
| Cryptographic event chain | WARNING: Planned | P2 | 11_event_logging.md |
| Audit verification tool | WARNING: Planned | P1 | 11_event_logging.md |
| Performance | | | |
| Database query threading | WARNING: Planned | P2 | 10_concurrency_model.md |

Tier 2 Summary: 13 features planned WARNING: (0 implemented)

Estimated Effort: 3-4 months (Phase 6)

TIER 3: Clinical Trials ([FAILED] NOT STARTED - Phase 7)

What This Enables: - Clinical trials with human subjects (IRB approval) - Real patient data collection - Treatment outcome tracking - Adverse event reporting

Requirements:

| Feature | Status | Priority | Notes |
|----------------------------------|----------------------|----------|-----------------------------|
| User Authentication | | | |
| Login system (username/password) | [FAILED] Not started | P0 | 08_security_architecture.md |
| Password hashing (Argon2) | [FAILED] Not started | P0 | 08_security_architecture.md |
| Session management | [FAILED] Not started | P0 | - |

| | | | |
|--|----------------------|----|-----------------------------|
| Password policy enforcement | [FAILED] Not started | P1 | 08_security_architecture.md |
| Access Control | | | |
| Role-based permissions | [FAILED] Not started | P0 | 08_security_architecture.md |
| Operator vs Administrator roles | [FAILED] Not started | P0 | 08_security_architecture.md |
| Audit trail for user actions | [FAILED] Not started | P0 | Already implemented [DONE] |
| Compliance | | | |
| HIPAA compliance validation | [FAILED] Not started | P0 | External audit |
| IRB documentation package | [FAILED] Not started | P0 | Clinical team |
| Informed consent tracking | [FAILED] Not started | P0 | New feature |
| Adverse event reporting | [FAILED] Not started | P0 | New feature |
| Data Protection | | | |
| PHI anonymization tools | [FAILED] Not started | P1 | New feature |
| Data export for analysis | [FAILED] Not started | P1 | New feature |
| Patient data retention policy | [FAILED] Not started | P0 | Policy document |
| Tier 3 Summary: 12 features required [FAILED] (0 started) | | | |
| Estimated Effort: 4-6 months (Phase 7) | | | |

TIER 4: FDA Submission/Production ([FAILED] NOT STARTED - Phase 8+)

What This Enables: - FDA 510(k) or PMA submission - Commercial production deployment - Clinical use outside of trials - Regulatory compliance across all jurisdictions

Requirements:

| Feature | Status | Priority | Notes |
|--|----------------------|----------|-------------------------|
| Design Controls | | | |
| Complete Design History File (DHF) | [FAILED] Not started | P0 | FDA requirement |
| Risk management file (ISO 14971) | [FAILED] Not started | P0 | FDA requirement |
| Traceability matrix (requirements ↔ tests) | [FAILED] Not started | P0 | 09_test_architecture.md |
| Design verification report | [FAILED] Not started | P0 | Test results summary |
| Design validation report | [FAILED] Not started | P0 | Clinical validation |
| Software Validation | | | |
| Software validation plan | [FAILED] Not started | P0 | IEC 62304 |
| Unit test documentation | [DONE] Implemented | - | 09_test_architecture.md |
| Integration test documentation | WARNING: Partial | P0 | Needs formal docs |
| System test documentation | [FAILED] Not started | P0 | End-to-end testing |

| | | | |
|--|----------------------|----|--------------------|
| Performance qualification (PQ) | [FAILED] Not started | P0 | Bench testing |
| Manufacturing | | | |
| Installation qualification (IQ) | [FAILED] Not started | P0 | Production units |
| Operational qualification (OQ) | [FAILED] Not started | P0 | Production units |
| Software bill of materials (SBOM) | [FAILED] Not started | P0 | Cybersecurity |
| Cybersecurity risk assessment | [FAILED] Not started | P0 | FDA requirement |
| Regulatory Submissions | | | |
| 510(k) pre-submission meeting | [FAILED] Not started | P0 | FDA interaction |
| 510(k) submission package | [FAILED] Not started | P0 | Complete filing |
| Post-market surveillance plan | [FAILED] Not started | P0 | FDA requirement |
| Quality System | | | |
| Quality manual (ISO 13485) | [FAILED] Not started | P0 | QMS documentation |
| Standard operating procedures (SOPs) | [FAILED] Not started | P0 | All processes |
| Change control procedures | [FAILED] Not started | P0 | Version management |
| CAPA system (Corrective/Preventive Action) | [FAILED] Not started | P0 | Quality system |

Tier 4 Summary: 21 features required [FAILED] (1 partially implemented)

Estimated Effort: 12-18 months (Phase 8+)

Quick Reference: What Can I Do?

[DONE] Right Now (Tier 1 - Current)

YES: - Develop and test algorithms - Run unit and integration tests - Benchmark performance - Train team members - Refine protocols - Test with calibration targets - Use mock patient data - Software development

NO: - Clinical trials - Real patient data - Animal testing - FDA submission - Production deployment

WARNING: After Phase 6 (Tier 2 - Pre-Clinical)

Additional YES: - Cadaver testing - Bench validation - Calibration verification - Safety system qualification - Animal testing (if applicable)

Still NO: - Human clinical trials - FDA submission - Production deployment

After Phase 7 (Tier 3 - Clinical Trials)

Additional YES: - IRB-approved human trials - Real patient data collection - Outcome tracking - Adverse event reporting

Still NO: - FDA submission (need full validation) - Commercial deployment

After Phase 8+ (Tier 4 - FDA/Production)

Additional YES: - FDA 510(k) submission - Commercial deployment - Clinical use - Production units

Development Roadmap

| CURRENT (Tier 1) | Phase 6 (Tier 2) | Phase 7 (Tier 3) | Phase 8+ (Tier 4) |
|--------------------------|---|--|---|
| Lab Testing 23 [DONE] | Pre-Clinical 13 WARNING: +3-4 months Encryption Calibration Audit Sigs | Clinical 12 [FAILED] +4-6 months Authentication Access Control IRB Compliance | FDA/Production 21 x +12-18 months Design Controls Validation Docs QMS Complete |

Total Features: - [DONE] 27 implemented (Tier 1 complete - updated Oct 2025) - WARNING: 13 planned (Tier 2 - Phase 6) - [FAILED] 33 future (Tiers 3 & 4 - Phases 7-8+) - 73 total features across all tiers

Risk Assessment by Tier

| Tier | Risk if Used Prematurely | Mitigation |
|---------------------|---------------------------------|---------------------------------|
| Tier 1 (Current) | Low - test data only | [DONE] No real patient exposure |
| Tier 2 (Phase 6) | Medium - unvalidated for humans | WARNING: Cadavers/bench only |
| Tier 3 (Phase 7) | High - regulatory violation | [FAILED] Need IRB approval |
| Tier 4 (Phase 8+) | Critical - illegal use | [FAILED] Need FDA clearance |

CRITICAL WARNING: - Using Tier 1 system for Tier 3 work (clinical trials) = HIPAA violation + IRB violation + potential FDA enforcement action - Using Tier 1 system for Tier 4 work (production) = Federal crime (selling unapproved medical device)

Summary

Current Capabilities (Tier 1): - Fully functional for lab experimentation and R&D - 27 core features implemented and tested (updated Oct 2025) - Safe for development with test data only - Production-ready camera implementation with Allied Vision API compliance

Next Milestone (Tier 2 - Phase 6): - 13 features to add (encryption, calibration, audit integrity) - Estimated 3-4 months development time - Enables pre-clinical validation and bench testing

Path to Clinical Use: - Tier 1 → Tier 2 → Tier 3 → Tier 4 - Total estimated time: 19-28 months from current state - Total features to add: 46 features across 3 phases

Document Owner: Project Manager Last Updated: 2025-10-30 Next Review: Start of Phase 6

Recent Updates (2025-10-30): - Added 4 new camera features (Allied Vision API integration, thread safety, hardware feedback, pixel format auto-detection) - Updated feature count from 23 → 27 implemented features - Updated FPS specification from 60 → 30 FPS (hardware-controlled frame rate) - Version: 1.1